

Lecture Notes: Glucose Sensor

1 Introduction

With over **346 million people** living with diabetes worldwide, glucose monitoring has become a vital part of modern healthcare.

Poor glucose management can lead to serious complications, and accurate detection of blood glucose is essential for diagnosis and control.

A **glucose sensor** is a **biosensor** that measures the glucose concentration in blood or biological fluids.

It consists of:

1. A **bioreceptor** (enzyme such as *glucose oxidase*).
2. A **transducer** that converts the biochemical signal into an electrical or optical one.
3. A **signal processor** that displays the glucose concentration digitally.

These sensors help in **continuous glucose monitoring** and are used in **glucometers**, **wearable devices**, and **clinical analyzers**.

2 Types of Glucose Sensors

According to your syllabus and textbook, the two major techniques used to examine glucose in blood are:

1. **Reflectance Photometry**
2. **Electrochemical Technique (Amperometric Biosensor)**

Note:

Although there is another electrochemical approach known as **Potentiometric sensing**, it is **not included as a main glucose sensor** in textbook because its **response depends strongly on pH and ion stability**, making it less reliable for practical glucose monitoring compared to reflectance photometry and amperometric methods.

Reflectance Photometric Glucose Sensor

Introduction

Reflectance photometry was one of the **earliest techniques** used in portable blood-glucose monitors.

It works on the principle of **light reflectance**, where the color intensity formed on a test strip after glucose reaction is measured.

Principle

- The **glucose oxidase enzyme (GOx)** on the test strip reacts with glucose in the blood sample.
 - The reaction produces a **colored complex** by oxidizing glucose into a molecule that reacts with a dye.
 - When light from an **LED** falls on this complex, the reflected light intensity changes based on glucose concentration.
 - A **photodiode** measures this reflected light, and a **microprocessor** converts it into a numerical glucose value.
-

Construction

- **Test strip:** coated with glucose oxidase and a **chromogenic indicator dye**.
 - **LED source:** emits light onto the strip.
 - **Photodiode:** detects reflected light.
 - **Microprocessor:** amplifies and interprets the signal to display the glucose level (in mmol/L).
-

Working

1. The patient pricks their finger using a **lancet device** to obtain a drop of blood.
2. The blood is applied to the enzyme-coated reagent pad.
3. Glucose in the blood reacts with glucose oxidase, forming a **colored compound**.
4. LED light reflects off the pad and reaches the **photodiode detector**.

5. The reflected intensity is converted into an **electrical signal**, processed by a **microprocessor**, and displayed as the glucose concentration.




Fig. Blood glucose monitor – Light from the LED reflecting from the reagent pad onto the photodiode

✓ Advantages

- Simple and convenient for quick tests.
- Low-cost and portable design.
- Fast visual measurement using disposable test strips.

⚠ Limitations

- Requires a **large blood volume** (1–3 μL).
- **Calibration** of test strips is essential.
- **Environmental factors** (light, humidity, and temperature) can affect readings.

4 Amperometric Glucose Biosensor

⚡ Introduction

Amperometric biosensors are the **modern electrochemical glucose sensors** found in most glucometers today.

They work by measuring the **electrical current** generated when glucose undergoes an


enzyme-catalyzed oxidation reaction.

This current is directly proportional to the glucose concentration in the sample.

Principle

The principle is based on **amperometry**, which measures the **current produced by redox reactions** at a fixed electrode potential.

The enzyme **glucose oxidase (GOx)** catalyzes the following reaction:



- The hydrogen peroxide (H_2O_2) produced is **oxidized electrochemically** at the electrode surface:
$$\begin{aligned} &[\\ &\text{H}_2\text{O}_2 \rightarrow 2\text{H}^+ + \text{O}_2 + 2\text{e}^- \\ &] \end{aligned}$$
- The electrons released during this oxidation flow through the external circuit, producing a **current (I)** that is proportional to the glucose concentration.

Thus, higher glucose levels produce a **stronger electrical current**.

Construction

Each **amperometric glucose sensor strip** contains multiple layers arranged in a compact structure:

1. **Working Electrode (Anode):**
 - Made of platinum, carbon, or gold.
 - Coated with immobilized *glucose oxidase enzyme*.
2. **Reference Electrode:**


- Usually **Ag/AgCl**, maintains constant potential.
3. **Counter Electrode:**
 - Completes the electrical circuit.
 4. **Spacer and Membrane Layers:**
 - Control glucose diffusion and block interfering chemicals.
 5. **Electronic Circuit and Display:**
 - Amplify, process, and display the current as glucose level.
-

Working

1. **Sample Application:** A small blood sample ($\sim 1 \mu\text{L}$) is applied on the test strip.
 2. **Enzymatic Reaction:** Glucose diffuses to the enzyme layer and is oxidized to gluconic acid and H_2O_2 .
 3. **Electrochemical Step:** The generated H_2O_2 reaches the electrode surface, where it is oxidized, releasing electrons.
 4. **Current Measurement:** The electrons travel through the circuit, generating a measurable current proportional to glucose concentration.
 5. **Signal Processing:** The device converts this current into a digital reading of glucose concentration, displayed in mmol/L or mg/dL.
-

Glucose Amperometric Biosensor (Glucometer)

Glucose oxidase to convert glucose to H_2O_2 and oxidizes it



(Use the "Glucose Amperometric Biosensor (Glucometer)" image from your book.)

✓ Advantages

- **Highly efficient** – gives results within 5 seconds.
- Requires **very small blood sample** (<1 μL).
- **More accurate and sensitive** than reflectance method.
- Suitable for **continuous monitoring** (used in CGM devices).

⚠ Limitations

- **Enzyme activity** decreases with time or high temperature.
- **Oxygen dependency** affects measurement in low- O_2 conditions.
- Requires **initial calibration** for consistent accuracy.

(✓ Matches your textbook perfectly — small blood volume, faster results, enzyme-based current generation.)

5 Comparison: Cost, Efficiency, Use, and Suitability

Aspect	Reflectance Photometric Sensor	Amperometric Sensor
Detection Type	Optical (light reflectance)	Electrochemical (current)
Principle	Measures light reflected from color formed by enzyme reaction	Measures oxidation current generated during glucose oxidation
Enzyme Used	Glucose Oxidase	Glucose Oxidase
Sample Volume	1–3 μL	< 1 μL
Response Time	15–20 seconds	4–6 seconds
Accuracy	Moderate	High
Cost	Low	Moderate–High
Where Used	Simple home test strips, educational kits	Digital glucometers and CGM systems
Why Used	Quick, low-cost screening	Precise, real-time diabetic monitoring
Main Limitation	Sensitive to external light and temperature	Sensitive to enzyme stability and oxygen
Overall Suitability	✓ Basic and affordable	🏆 Best for accurate, fast, and practical use

6 Conclusion

- **Reflectance photometric sensors** introduced the basic optical principle for glucose testing — simple but calibration-dependent.
- **Amperometric biosensors** advanced glucose monitoring with high speed, sensitivity, and precision, becoming the standard in all modern glucometers.
- **Potentiometric sensors**, though scientifically valid, are not included as a main type in your syllabus because their **voltage response depends on pH and ion stability**, making them less practical for real-world glucose measurement.

💡 Summary Thought:

“Glucose sensors represent the perfect blend of chemistry, biology, and electronics — a life-saving application of engineering principles that bridges science with human health.”

