

Excellent — that's a very teacher-minded request 🙌

Here's a set of **likely student questions and clear, confident answers** based on your *Glucose Sensor* lecture.

They're written in a **Q&A format**, covering **conceptual**, **application**, and **clarification**-type doubts that students would realistically ask during or after your session.

You can use these for quick discussion, viva prep, or tutorial sheets.



Possible Student Doubts — Glucose Sensor Lecture



1 What is the main difference between a biosensor and a normal sensor?

Answer:

A **biosensor** uses a **biological element** (like an enzyme, antibody, or microorganism) to detect a specific substance and convert it into a measurable signal.

A **normal sensor** only detects physical quantities (like temperature or pressure) using physical transducers — it doesn't involve a biological reaction.



2 Why is Glucose Oxidase (GOx) used in glucose sensors?

Answer:

Glucose Oxidase specifically reacts **only with glucose molecules**, converting them into **gluconic acid and hydrogen peroxide (H_2O_2)**.

This selective nature makes it ideal for accurate glucose detection without interference from other sugars.



3 How does the reflectance photometric sensor actually detect glucose?

Answer:

It measures the **intensity of light reflected** from a **color-changing reagent pad** on the test strip.

When glucose reacts with the enzyme and dye, the color intensity changes — the **darker the color**, the **higher the glucose level**.

A **photodiode** senses this reflected light, and a **microprocessor** converts it to a glucose value.



4 Why does the color become darker when glucose concentration increases?

Answer:

More glucose means a stronger enzyme reaction, producing more of the **colored dye complex**.

This absorbs more light and reflects less — so the pad looks **darker**, and the sensor detects a lower reflectance signal (indicating high glucose).




5 In amperometric sensors, why does the current increase with glucose level?

Answer:

Because more glucose molecules produce more **hydrogen peroxide (H₂O₂)** during the enzymatic reaction.

Each H₂O₂ molecule releases **two electrons** when oxidized at the electrode:



The greater the number of electrons, the **higher the current**, so current is **directly proportional** to glucose concentration.



6 What does “amperometric” actually mean?

Answer:

“Amperometric” comes from **ampere**, the unit of electric current.

It means that the sensor measures **current** produced by oxidation or reduction reactions — in this case, during the oxidation of glucose or hydrogen peroxide.



7 What is the purpose of the reference electrode in an amperometric sensor?

Answer:

The **reference electrode** (usually Ag/AgCl) maintains a **constant potential** so the working electrode can operate at a fixed voltage.

This ensures that the measured current depends **only on glucose concentration**, not on voltage fluctuations.



8 What happens if the enzyme Glucose Oxidase becomes inactive?

Answer:

If the enzyme is denatured (by heat, aging, or chemicals), it **can't catalyze glucose oxidation**, so the sensor produces **no reliable signal**.

That's why enzyme-based sensors have an **expiry date** and must be stored properly.



9 Why do we need membranes in amperometric biosensors?

Answer:

Membranes act as **selective barriers** — they allow **glucose** to diffuse in but block **interfering compounds** (like ascorbic acid, uric acid).

They also **control the diffusion rate** to ensure stable readings and prevent enzyme leakage.



10 Why is hydrogen peroxide (H_2O_2) important in the amperometric method?

Answer:

Hydrogen peroxide is the **electro-active product** formed during glucose oxidation.

Its oxidation at the electrode releases electrons, which generate the **measurable current**.

Without H_2O_2 formation, there would be no current flow for measurement.



1 1 How does the sensor convert current into glucose concentration numerically?

Answer:

The microprocessor uses **calibration data** (pre-set correlation between current and glucose concentration).

When current is measured, it compares it with stored calibration curves and displays the corresponding glucose value in **mmol/L** or **mg/dL**.



1 2 Why are potentiometric glucose sensors not used widely?

Answer:

Because their **voltage (potential) output** depends on **pH and ion stability**.

Small variations in blood pH can cause large voltage changes, making them **less accurate and less stable** compared to amperometric sensors.



1 3 Which sensor type is faster — Reflectance or Amperometric? Why?

Answer:

The **Amperometric sensor** is faster (4–6 seconds) because the **electron transfer reaction** occurs instantly after glucose oxidation.

Reflectance sensors require **color development time**, usually 15–20 seconds.



14 Why does the amperometric sensor require less blood than the photometric one?

Answer:


Because the electrochemical reaction needs only a **thin layer of sample** to complete the redox process, whereas the photometric method needs a **larger visible area** for accurate color detection.



15 What is the role of oxygen in the enzyme reaction?

Answer:

Oxygen acts as the **electron acceptor** in the oxidation of glucose:



It helps in producing **hydrogen peroxide**, which is later oxidized at the electrode.



16 Why do some advanced biosensors use mediators like ferrocene?

Answer:

Mediators are **electron carriers** that transfer electrons between the enzyme and electrode more efficiently than oxygen.

This allows sensors to work **even at low oxygen levels** and improves sensitivity — these are called **second-generation amperometric sensors**.



17 What units are used to express glucose concentration?

Answer:

Usually in **mg/dL (milligrams per deciliter)** or **mmol/L (millimoles per liter)**.

Conversion:

$1 \text{ , } \text{mmol/L} = 18 \text{ , } \text{mg/dL}$



18 Why is calibration important in both sensors?

Answer:

Calibration ensures the **measured signal (light or current)** correctly matches known glucose concentrations.

Without calibration, the displayed values can **drift or be inaccurate** due to aging of components or environmental changes.



19 Can glucose sensors measure other sugars?

Answer:

Not effectively — **Glucose Oxidase (GOx)** is specific only to **β -D-glucose**.

Other sugars (like fructose or galactose) do not react significantly, so they have negligible effect on readings.



20 What are the future developments in glucose sensing?

Answer:

- **Non-invasive sensors** using **infrared spectroscopy** or **fluorescent methods**.
- **Nanomaterial-based electrodes** for higher sensitivity.
- **Enzyme-free sensors** using catalytic metals like platinum or graphene.
- **Continuous monitoring patches** integrated with mobile devices.



Quick Summary of Common Viva Questions

Concept	Typical Viva Question	Short Answer
Enzyme used	Which enzyme is used in glucose sensors?	Glucose Oxidase (GOx)
Basic Principle	What is the principle of glucose sensing?	Enzymatic oxidation of glucose measured as optical or electrical signal
Product of Reaction	What are the products formed in glucose oxidation?	Gluconic acid and hydrogen peroxide (H_2O_2)
Sensor Type	What are the two main glucose sensors?	Reflectance photometric and amperometric biosensors
Electrochemical Basis	What does amperometric mean?	Measuring current due to redox reaction
Potentiometric Exclusion	Why not potentiometric?	Voltage depends on pH and ion interference
Time Difference	Which is faster?	Amperometric (~5 s)
Sample Volume	Which needs more blood?	Reflectance photometric (~1–3 μ L)

Concept	Typical Viva Question	Short Answer
Application	Where is amperometric used?	Clinical and wearable glucose monitoring

✅ **Tip for Teaching or Viva:**

If students ask, “Which sensor is better?”, you can conclude confidently:

“Reflectance photometric sensors were the starting point, but amperometric biosensors are the gold standard today — faster, smaller, more accurate, and ideal for continuous glucose monitoring.”

Would you like me to generate a **short 10-question multiple-choice quiz** (with answers and reasoning) that you can use to test students after your lecture or for self-practice?